

# Sol-Cart B<sup>®</sup>

## Product Information

### Cartridge with sodium bicarbonate for hemodialysis

#### Bicarbonate cartridge

- Cartridge housing material: polypropylene
- Filter: one filter on the inlet and one filter on the outlet port of the cartridge
- Volume of contents: 650 g / 760 g / 1100 g

#### Composition

- Sodium bicarbonate powder meeting the monograph of the European Pharmacopeia (medical device)

#### Compatibility

- Sol-Cart B<sup>®</sup> fits to all current dialysis machines with a suitable cartridge holder

#### Presentations

	Cartridges per box	Art-No	Cartridges per pallet
■ Sol-Cart B <sup>®</sup> 650 g	10 x 650 g	496	630
■ Sol-Cart B <sup>®</sup> 760 g	10 x 760 g	494	630
■ Sol-Cart B <sup>®</sup> 1100 g	8 x 1100 g	804	432



**B. Braun Avitum AG**

Schwarzenberger Weg 73–79  
34212 Melsungen

Germany

Tel +49 (56 61) 71-26 24

Fax +49 (56 61) 75-26 24

e-mail: [dialysis@bbraun.com](mailto:dialysis@bbraun.com)

[www.bbraun-dialysis.com](http://www.bbraun-dialysis.com)